LABIOSERVICE

Labioservice S.r.l.

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C.F. - 02951070586 / **P.IVA** - 01140091008

Technical Specifications

Name	Value
Packaging Contents	1x RF542 (including 3.6V AA Lithium Battery), 1x RF542 Wall Mount Bracket, 2x Cable Ties, 1x RF520 AC Adapter Including Regional Adapter, 1x Cat5e Ethernet Cable, 1x RF525 Transmitter Activator, 1x RF542 Quick Reference Guide
Operating Temperature Range	0°C to 50°C 10-90%RH Non-condensing
Temperature Measurement Range - External Sensor	-40°C to +125°C / - 40°F to +257°F
System Accuracy with External Thermistor - 20°C to +70°C	±0.5°C / ±0.9°F
System Accuracy with External Thermistor - Full Range	±1°C / ±2°F
Resolution Temperature - Thermistor	0.1°C / 0.2°F
Door Sensor	7.5 seconds
Storage Temperature	-40°C to +85°C / -40°F to +185°F
Ethernet Connectivity	Fixed 100Mbps Half Duplex RJ45
Network Configuration	DHCP
Clock Accuracy	20ppm (1 minute/month) at 25°C / 77°F
Logging Memory	32,000 records
Logging Frequency	Programmable between 1 minute and 60 minutes
Monitoring Frequency	1 minute
LEDs	Red (Front Label) - Alarm Green (Front Label) - External Power Green (Side) - LAN Connection
Case Material	Over Moulded food safe clear Polycarbonate with BioCote® antimicrobial
Battery Type	Replaceable Lithium AA Cell 3.6V
Battery Life	Up to 2 years backup, Normal operation requires RF520 (Supplied)
Dimensions	L 83mm x W 107mm D 40mm
Weight	215g

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
C.F. - 02951070586 / **P.IVA -** 01140091008

Name	Value
Probe Attachment Required	A standard Comark probe is required for operation of this device
Warranty	2 Years

Battery Change Process

To change the battery in your RF542 Remote Temperature Monitor simply follow the following steps:

- Remove mains power from the RF542 Remote Temperature Monitor
- Undo the battery cover and remove the Lithium AA battery
- Wait for the screen to go off completely
- Insert the new RF542BATT and re-fit the battery cover
- Re-connect the mains power
- The RF542 Remote Temperature Monitor should re-connect to the Gateway automatically
- Dispose of the old battery safely – Note Lithium Batteries should not be disposed of in General Waste

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Technical Specifications

Battery Type AA, Lithium Thionyl Chloride Battery (Li-SOCI2)
Nominal Voltage 3.6 V
Nominal Capacity 2.2 Ah